

# Roche SARS-CoV-2 Rapid Antigen Test (Nasopharyngeal swab)

In times of uncertainty and increasing patient numbers, the reliable SARS-CoV-2 Rapid Antigen Test provides fast answers, wherever you need them.



## Specifications

- Efficient and cost-friendly
- Results in 15 minutes
- Storage Condition: 2-30°C/36-86°F
- Sensitivity: 86.2 - 90.0%
- Specificity: 97.9 - 99.3%
- Intended for use with Nasal and Nasopharyngeal swab specimens within 5 days of symptom onset



## Packaging

- Quantity: 25 Test Devices (Packed Individually)
- Instructions for Use
- Quick Reference Guide
- Extraction buffer tube
- Nozzle cap
- Sterile Swab
- Film