

### **Product Description**

Product name: KN95-04 protective mask.

Main structure: Composed of non-woven fabric with folded ultrasonic compound mask body, nose clip and mask band.

Product Size: 150x100mm(10pcs/box)

Instructions for use:

Put the colored side of the mask outwards, Cover nose and mouth, Hang the loop on the ear to adjust the radian of the film on the

nose bridge to fit snugly,

Pull the mask down to cover your chin.

#### Precautions:

1. This product is for single use only and it is forbidden to reuse it.

2.After use, do not discard it randomly, you should discard a maskspecific trash can.

#### Storage conditions:

Should be stored in a cool, dry place and keep away from heat sources, avoid direct sunlight.

Guideline: GB2626-2006

**Expiration date:** Stored under specified conditions, valid for 24

months.

10pcs/box







### 3D Stereo Tailoring Design

All-round fit to facial contours unbound effectively isolate pollution













## Quadruple Structure Security Protection

Quadruple structure high-efficiency electrostatic filter protects the respiratory tract not affected by yellow sand and dust.









# Mask Packaging





Gift box size: 180X120X40mm (10 pcs/box)

G.W.: 0.091Kg

Master Carton size: 510X380X440mm (80 pack/CTN)

G.W.: 8.6Kg

Product real shot display











#### **Certificates**



